

## **E. coli growth media**

### **LB:**

For 1 L

- 10 g tryptone
- 5 g yeast extract
- 10 g NaCl
- Optional: Bring up the volume to around 900 mL with ddH<sub>2</sub>O, then adjust the pH to 7.4
- Bring up the volume to 1 L with ddH<sub>2</sub>O
- Sterilize by autoclaving for 30 min

To make LB-agar, add 15 g of agar prior to autoclaving

### **Low-salt LB:**

For 1 L

- 10 g tryptone
- 5 g yeast extract
- 5 g NaCl
- Optional: Bring up the volume to around 900 mL with ddH<sub>2</sub>O, then adjust the pH to 7.4
- Bring up the volume to 1 L with ddH<sub>2</sub>O
- Sterilize by autoclaving for 30 min

To make low-salt LB-agar, add 15 g of agar prior to autoclaving

### **SOC:**

For 1 L

- 20 g tryptone
- 5 g yeast extract
- 0.5 g NaCl
- 0.186 g KCl
- 0.952 g MgCl<sub>2</sub>
- Bring up the volume to around 900 mL with ddH<sub>2</sub>O, then adjust the pH to 7.4 with 10 *N* NaOH
- Bring up the volume to 1 L with ddH<sub>2</sub>O
- Sterilize by autoclaving for 30 min
- Add 20 mL of sterile 1 M glucose immediately before use

### **1 M glucose:**

For 100 mL

- Dissolve 18 g glucose in 100 mL ddH<sub>2</sub>O
- Sterilize by filtering through 0.2 μ or by autoclaving for 15 min

## **Antibiotic stocks**

### **Ampicillin (1000x):**

- Dissolve 5 g ampicillin in 25 mL ddH<sub>2</sub>O
- Add 25 mL absolute ethanol
- Store at -20 °C

**Carbenicillin (1000x for liquid media, 2000x for plates):**

- Dissolve 2.5 g carbenicillin in 25 mL ddH<sub>2</sub>O
- Add 25 mL absolute ethanol
- Store at -20 °C

**Chloramphenicol (1000x):**

- Dissolve 1.7 g chloramphenicol in 50 mL absolute ethanol
- Store at -20 °C

**Kanamycin (1000x, but 250x for autoinduction media):**

- Dissolve 1.25 g ampicillin in 50 mL ddH<sub>2</sub>O
- Sterilize by filtering through 0.2 μ
- Store at -20 °C in 1.5-mL aliquots

**Tetracycline (1000x):**

- Dissolve 0.75 g tetracycline in 50 mL absolute ethanol or isopropanol
- Store at -20 °C in 1.5-mL aliquots