

## Transformation of *E. coli*

### Use aseptic techniques!

#### Calcium chloride method:

For pure plasmid DNA clones, transform 10 ng of plasmid into 10  $\mu\text{L}$  of competent cells. For ligation reactions, use 10  $\mu\text{L}$  of cells for each  $\mu\text{L}$  of ligation mix.

- (1) Thaw the competent cells on ice.
- (2) Mix the DNA and competent cells in a sterile 1.5-mL microfuge tube.
- (3) Incubate in ice for 30 min.
- (4) Heat shock at 42 °C for 30 sec (1 min for volumes  $\geq 50 \mu\text{L}$ ).
- (5) Incubate in ice for 2 min.
- (6) Add 9 volumes of SOC or LB medium containing **no antibiotics**. Media pre-warmed to 37 °C work best, but this is not required. Shake for at least 30 min (45-60 min if the plasmid encodes resistance to an antibiotic other than ampicillin). When transforming pure plasmids encoding ampicillin resistance, steps 6 and 7 are not required.
- (7) If the total volume is  $>100 \mu\text{L}$ , spin down the cells in a microcentrifuge for 1 min at 8,000 g. Remove all but  $\sim 100 \mu\text{L}$  of the supernatant. Resuspend the cells in the remaining liquid.
- (8) Plate on LB-agar (preferably pre-warmed) containing the appropriate antibiotics.
- (9) Allow the plate to dry for around 10 mins, then incubate upside-down at 37 °C overnight.

#### Things to look up / think about:

How does calcium chloride make *E. coli* competent to take up DNA?

Why do we pre-grow the cells in media with no antibiotics prior to plating? Why can we skip this step if the antibiotic is ampicillin?